




### Determination of Caffeine in Soda Using HALO® C18, 5 µm

Application Note 145-F



Sample	Caffeine tested mg/(355 mL)	Can value mg/(355 mL)
Store brand cola 1	12	N/A
Cola 2	53	54
Cola 3	43	43
Cola 4	36	38
Cola 5	38	38
Store brand diet cola 1	12	N/A
Diet cola 2	45	46
Diet cola 3	34	34
Diet cola 4	36	35
Energy drink 1*	160	160
Energy drink 2**	79	80
Diet Energy drink**	79	80
Non-cola drink 1	53.3	54
Non-cola drink 2	22	22
Diet non-cola drink	43	41
Diet cola 1 non-caffeinated	0	N/A
Diet cola 2 non-caffeinated	0	N/A
Diet cola 3 non-caffeinated	0	N/A

355 mL = 12 oz.

\*amount in 16 oz. (473 mL) cans

\*\*amount in 8.4 oz (248 mL) cans

Caffeine is a stimulant found at various levels in coffee, colas, and energy drinks. HPLC is a convenient way to determine the amount of caffeine present. Here, sodas were analyzed by direct injection onto a 5 µm HALO® C18 column after decarbonation. A guard column should be used in this application.

#### TEST CONDITIONS:

**Column:** HALO 90 Å C18, 5 µm, 3.0 x 50 mm,  
HALO 5 µm guard column

**Part Numbers:** 95813-402, 95813-102

**Mobile Phase:** 75/25 - A/B

A: 0.1% formic acid in water

B: Methanol

**Flow Rate:** 0.8 mL/min

**Pressure:** 120 bar

**Temperature:** 30 °C

**Detection:** UV 254 nm, VWD

**Injection Volume:** 1.0 µL

**Sample Solvent:** (Caffeine std.) mobile phase

**Response Time:** 0.02 sec


**Data Rate:** 25 Hz

**Flow Cell:** 2.5 µL semi-micro

**LC System:** Shimadzu Prominence UFLC XR

**Extra Column Volume:** ~14 µL

#### STRUCTURE:



Caffeine

