




## Separation of Neutral Aromatics on HALO® PFP, C18 and Phenyl-Hexyl

Application Note 23-N



### PEAK IDENTITIES:

1. Butylbenzene
  2. Acenaphthene
  3. 1-Phenyl-naphthalene
  4. Pyrene
- i = impurities

The separation of nonpolar aromatic compounds on these three HALO® bonded phases under the same conditions show differences in selectivity that can be utilized in optimizing difficult separations.

### TEST CONDITIONS:

#### Columns:

- 1) HALO 90 Å PFP, 2.7  $\mu\text{m}$ , 4.6 x 50 mm  
Part Number: 92814-409
- 2) HALO 90 Å C18, 2.7  $\mu\text{m}$ , 4.6 x 50 mm  
Part Number: 92814-402
- 3) HALO 90 Å Phenyl-Hexyl, 2.7  $\mu\text{m}$ , 4.6 x 50 mm  
Part Number: 92814-406

#### Mobile Phase: 30/70 - A/B

A: Water

B: Methanol

**Flow Rate:** 2.0 mL/min

**Pressure:** ~250 bar

**Temperature:** 40 °C

**Detection:** UV 254 nm, VWD

**Injection Volume:** 1.0  $\mu\text{L}$

**Sample Solvent:** Methanol


**Response Time:** 0.02 sec

**Flow Cell:** 2.5  $\mu\text{L}$  semi-micro

**LC System:** Shimadzu Prominence UFLC XR

**Extra Column Volume:** ~14  $\mu\text{L}$


### STRUCTURES:




Butylbenzene



1-Phenyl-naphthalene



Acenaphthene



Pyrene

