




## High Temperature/Low pH Stability with HALO 160 Å ES-C18, 2.0 µm

Application Note 137-PE



The sterically-protected C18 phase on the HALO® 2.0 µm 160 Å column enables high temperature stability with low pH mobile phases. The replicate injections were stopped at injection 480 (15,500 column volumes). The column is expected to have a lifetime of ~1000 injections, depending on the type of sample and conditions used.

### PEAK IDENTITIES: MW (g/mol):

1. Gly-Tyr	238
2. Val-Tyr-Val	380
3. Met-enkephalin	574
4. Angiotensin II	1046
5. Leu-enkephalin	556
6. Ribonuclease A	13,700
7. Bovine insulin	5733

### TEST CONDITIONS:

**Column:** HALO 160 Å ES-C18, 2.0 µm,  
2.1 x 100 mm  
**Part Number:** 91122-602  
**Mobile Phase:**  
A: 0.1% trifluoroacetic acid in water  
B: 0.1% trifluoroacetic acid in 80/20 acetonitrile/  
water  
**Gradient:** 6% B to 54% B in 10 min  
**Flow Rate:** 0.5 mL/min  
**Initial Pressure:** 395 bar  
**Maximum Pressure:** 417 bar  
**Temperature:** 60 °C  
**Detection:** UV 215 nm, PDA  
**Injection Volume:** 0.5 µL  
**Sample Solvent:** Mobile phase A  
**Response Time:** 0.025 sec  
**Data Rate:** 40 Hz  
**Flow Cell:** 1.0 µL  
**LC System:** Shimadzu Nexera X2

