



New Reversed-Phase SPP Columns with Alternate Selectivity for Small Molecules

Richard A. Henry, Independent Consultant, State College, PA 16801

Stephanie A. Schuster, Conner W. McHale,

William L. Johnson, Joseph J. DeStefano

Advanced Materials Technology, Inc., Wilmington, DE 19810


Outline

1. Description of Fused-Core[®] SPP silica for small molecule HPLC.
2. HALO[®] AQ-C18: designed with adequate polarity for pores to stay wet in pure water; selectivity toward highly polar analytes is altered very slightly.
3. HALO[®] Biphenyl: designed with increased aromatic character to improve retention and selectivity for polar analytes over less polar phases like Phenyl-Hexyl.
4. Discussion of column selectivity concepts, F_s equivalency values, and the USP/PQRI database.
5. Selected applications to illustrate and compare selectivity.

The Pioneering Shell Particle for Small Molecules


Fused-Core[®] Specifications:

- Superficially Porous (SPP Type)
- Type B silica (highly pure and low acidity)
- Available in 2, 2.7, and 5 μm total size*
- Porous shell with controlled thickness for short diffusion paths
- 90 \AA pores for unrestricted diffusion and high surface area for good retention and high loadability of small molecules




* 2 and 5 μm Biphenyl available in 2019

New HALO[®] AQ-C18 and Biphenyl Phases



AQ-C18



Biphenyl

- Designed to allow phase pores to stay wet in low organic and to add a small, reproducible polar selectivity component
- Designed to change aromatic/aliphatic content and alter selectivity for many polar analytes

F_s Values for 90 Å HALO[®] Phases vs HALO[®] C18

Columns shown in descending order of hydrophobicity, H.

F_s	Column	H	S^*	A	B	C (pH 2.8)	C (pH 7.0)	ethylbenzene k
0	HALO C18	1.1	0.04	0	-0.05	0.05	0.04	6.1
12.07	HALO AQ-C18	1	-0.036	0.099	-0.048	0.156	0.864	6.7
10.04	HALO C8	0.91	0.02	-0.13	0	-0.01	0.18	4.3
52.83	HALO RP-Amide	0.85	0.08	-0.38	0.19	-0.41	0.31	4.6
17.35	HALO Phenyl-Hexyl	0.78	-0.09	-0.23	0	0.1	0.45	3.5
26.76	HALO Biphenyl	0.708	-0.183	-0.279	0.028	0.047	0.99	3.1
94.45	HALO PFP	0.702	-0.117	-0.073	-0.062	1.17	0.972	2.3
22.78	HALO ES-CN	0.566	-0.11	-0.344	0.021	0.126	1.15	1.88

$$F_s = \sqrt{(w_H(\mathbf{H}_1 - \mathbf{H}_2))^2 + (w_{S^*}(\mathbf{S}_1^* - \mathbf{S}_2^*))^2 + (w_A(\mathbf{A}_1 - \mathbf{A}_2))^2 + (w_B(\mathbf{B}_1 - \mathbf{B}_2))^2 + (w_{C_{2.8}}(\mathbf{C}_{2.81} - \mathbf{C}_{2.82}))^2}$$


Data provided online by Stoll group at Gustavus Adolphus College. www.hplccolumns.org

Overview of Two New HALO[®] HPLC Columns

- HALO[®] AQ-C18
 - Maintains L1 column classification (>90% C18 character).
 - Complements HALO[®] C18 (in ease of use, stability and reproducibility).
 - Moderate change in retention and selectivity for many polar analytes.
 - Useful in both methanol and acetonitrile mobile phases.
 - Extends mobile phase range to 100% aqueous and low organic for polar analytes, including gradients.
- HALO[®] Biphenyl
 - Maintains L11 USP column classification (phenyl aromatic phase).
 - Complements HALO[®] Phenyl-Hexyl with alternate selectivity.
 - Shows significant selectivity differences for polar and aromatic analytes when compared to phenyl phases with alkyl character.
 - Very useful in both methanol and acetonitrile mobile phases.


Purines on Aqueous-Compatible C18 SPP Phases

2.1 x 100 mm, 0.1% DFA in water, 0.35 mL/min, 35 °C, 1 µl @ 10 mM each in 0.1% DFA, UV detection @ 265 nm




2.7 µm HALO 90 Å AQ-C18

1. Guanine
2. Hypoxanthine
3. Uracil
4. Xanthine




2.7 µm SPP Aqueous Compatible Phase



2.6 µm SPP Aqueous Compatible Phase

Comparative results presented here may not be representative for all applications.

Comparison of Aromatic Test Mix in Acetonitrile



Columns: HALO 90 Å C18, 2.7 µm; HALO 90 Å AQ-C18, 2.7 µm; HALO 90 Å Biphenyl, 2.7 µm; all 4.6 x 100 mm

Mobile Phase A: water

Mobile Phase B: **ACN**

Gradient:	Time	%B
	0	40
	15	90

Flow Rate: 1.5 mL/min

Temperature: 30 ° C

Detection: UV 254 nm, VWD

Injection Volume: 2 µL

Peak Identities:

1. uracil
2. resorcinol
3. 4-nitrophenol
4. benzonitrile
5. anisole
6. trans-stilbene
7. valerophenone
8. biphenyl
9. impurity
10. o-terphenyl

Comparison of Aromatic Test Mix in Methanol

Columns: HALO 90 Å C18, 2.7 µm; HALO 90 Å AQ-C18, 2.7 µm; HALO 90 Å Biphenyl, 2.7 µm; all 4.6 x 100 mm

Mobile Phase A: water

Mobile Phase B: **methanol**


Gradient:	Time	%B
	0	45
	20	55

Flow Rate: 1.2 mL/min

Temperature: 30 ° C

Detection: UV 254 nm, VWD





Injection Volume: 2 µL



Peak Identities:


1. uracil
2. resorcinol
3. 4-nitrophenol
4. benzonitrile
5. anisole
6. valerophenone
7. trans-stilbene
8. biphenyl
9. o-terphenyl
10. impurity

Selectivity Comparison for Aromatic Test Mix



See previous two slides for conditions

Comparison of Sulfa Drugs on Four Phases



Columns: HALO 90 Å C18, 2.7 µm; HALO 90 Å AQ-C18, 2.7 µm; HALO 90 Å Phenyl-Hexyl, 2.7 µm; HALO 90 Å Biphenyl, 2.7 µm; all 3.0 x 100 mm

Mobile Phase A: 20 mM Ammonium Formate, pH 3

Mobile Phase B: **methanol**

Gradient:	Time	%B
	0	12
	1.5	13
	4.0	50
	6.0	60

Flow Rate: 0.65 mL/min

Temperature: 30 ° C

Detection: UV 260 nm, VWD

Injection Volume: 1 µL

Peak Identities:

1. uracil
2. sulfadiazine
3. sulfathiazole
4. sulfamerazine
5. sulfamethizole
6. sulfachloropyridazine
7. sulfamethoxazole
8. sulfadimethoxin

Estrogen Steroid Separation on Three Phases

Mobile Phase A: 10mM Potassium Phosphate pH: 7

B: Acetonitrile

Gradient: 20-60%B in 6 min.

Instrument: Nexera

Wavelength: 220nm

Injection: 1 μ l


Temperature: 25°C

Flow Rate: 0.4 mL/min.


Column: 2.1x100mm, 2.7 μ m

Peak Identities

1. B-Estradiol-3-sulfate
2. Estriol
3. Estrone-3-sulfate
4. B-estradiol
5. A-estradiol
6. Androstenedione
7. Estrone



Steroid Separation on HALO[®] Biphenyl



PEAK IDENTITIES:

1. Estriol
2. Hydrocortisone
3. Prednisone
4. Cortisone
5. Corticosterone
6. β -Estradiol
7. Cortisone Acetate
8. Testosterone
9. 17- α -Hydroxyprogesterone
10. 11-Deoxycorticosterone
11. Progesterone

Column: HALO 90Å Biphenyl, 2.7 μ m, 4.6 x 50mm

Mobile Phase A: Water

Mobile Phase B: Acetonitrile

Gradient: 20-60% B in 6 minutes

Flow Rate: 1.85 mL/min

Initial Pressure: 344 bar

Temperature: 30 °C


Detection: UV 215 nm, PDA

Injection Volume: 4 μ L

Sample Solvent: Acetonitrile:water, 37.5:62.5

LC-MS of Pesticides on Biphenyl

HALO® Biphenyl



Column: HALO 90 Å Biphenyl, 2.7 µm, 2.1 x 100 mm

Flow Rate: 0.2 mL/min

Injection Volume: 1 µL

Sample: Nevada Pesticide Mix

Detection: MS

Mobile Phase: A= Water/0.1% formic acid/4 mM ammonium formate

B= Acetonitrile/0.1% formic acid/4 mM ammonium formate

Gradient:	Time (min.)	%B
	0.00	0
	1.01	15
	4.00	35
	5.00	62
	30.00	100
	34.00	100

18,19

PEAK IDENTITIES

1. Daminozide
2. Flonicamid
3. Thiamethoxam
4. Imidacloprid
5. Pacllobutrazol
6. Fenhexamid
7. Myclobutanil
8. Bifenazate
9. Dimethomorph Isomer 1
10. Pirotetramat
11. Dimethomorph Isomer 2
12. Spinosad A
13. Spinosad D
14. Trifloxystrobin
15. Spinetoram
16. Pyrethrin II
17. Piperonyl butoxide
18. Pyrethrin I
19. Etoazole
20. Abamectin A
21. Cypermethrin
22. Bifenthrin
23. Acequinocyl

1

2

3 4

5

7

6

8

11

9

12

13,14

15,16,17

20 21

22

23

Summary and Conclusions

- HALO[®] AQ-C18 allows exploration of a wider solvent range starting with low organic or even 100% aqueous without concern for stationary phase dewetting; major change in C18 selectivity or retention is unlikely because the C18 character greatly dominates.
- Both MeCN and MeOH should be evaluated when HALO[®] AQ-C18 and Biphenyl columns are selected for method development; closely separated analytes with polar character may increase in resolution and in some cases may even change elution order in MeOH.
- HALO[®] Biphenyl with dominant aromatic character shows enhanced interactions with polar groups compared to either C18 or AQ-C18 and has even greater aromatic character than HALO[®] Phenyl-Hexyl, which has about equal aliphatic and aromatic content.

Acknowledgements

- Tim Langlois and Barry Boyes
- Ben Libert, Justin Godinho and Chuping Luo
- Matt Jackson and Brian Wagner

Trademarks:

Fused-Core and HALO are registered trademarks of Advanced Materials Technology, Inc.

www.advanced-materials-tech.com

www.mac-mod.com